

# What's new in AnywhereDolphin

## C O N T E N T S :


|  |   |
|--|---|
| NEW RESPONSIVE DESIGN .....  | 2 |
| REMOVED DEPENDENCY ON MICROSOFT<br>SILVERLIGHT PLUGIN .....                | 3 |
| THE FREE DOLPHIN IMAGE VIEWER NOW<br>INTEGRATES WITH ANYWHERE DOLPHIN..... | 4 |

---

## New Responsive Design

---

To respond to the increased access to AnywhereDolphin from mobile devices, we've re-engineered AnywhereDolphin to responsive design. Responsive simply means AnywhereDolphin sizes to the screen size it is being displayed. This improves the usability and experience on smaller screens, like phones and tablets.



---

## Removed Dependency on Microsoft Silverlight Plugin

---

There were a couple components of AnywhereDolphin, most notably, the tooth chart and Ceph Viewer, that relied on the Microsoft Silverlight plugin. Microsoft has discontinued future updates to Silverlight and it is recommended to remove future dependency to Silverlight. We've been hard at work the last couple of months taking care of this.

You'll not see any changes in functionality, but they look and feel will be different.




---

## The free Dolphin Image Viewer now integrates with AnywhereDolphin

---

You can now download images shared with AnywhereDolphin using the free Dolphin Viewer. This makes it easier to share records, including 3D files with non-Dolphin Imaging customers. From the viewer, simply enter your AnywhereDolphin credentials to access patient records shared via AnywhereDolphin:



You can read more about the Image Viewer and download it from here: The free viewer can be downloaded from here: <http://www.dolphinimaging.com/Product/ImagingViewer>