



Dolphin 3D

Import your digital study models into Dolphin!

The new Digital Study Model software module lets you store and organize your 3D digital study model files! Whether you scan your models or take a digital impression, you can use Dolphin to store these records within the patient's chart.

- Supports standard STL file formats
- 360° view of models within Dolphin with several layout options
- Patient CBCT is not required to store models
- Superimpose "teeth" onto a patient's CBCT scan for viewing and analysis*

* *Dolphin 3D is required to superimpose "teeth" onto patient's CBCT scan.*

